

INSTRUCTIONS

ALPHASEAL UNIVERSAL CAVITY LINER With Fluoride Light-Cured

Kit Contains
12 g of Liner
in Dentin Shade
Instructions

GENERAL INFORMATION

Alphaseal is a visible light cured cavity liner. It has a light paste consistency to facilitate the application. This material consists of a biocompatible resin, filled with hydroxyapatite, barium sulfate, and fluoride salts. Alphaseal is compatible with all commercially available restorative materials and cements. When used under water-based cements, however, wiping the surface of cured Alphaseal with a cotton pledget moistened with alcohol is recommended, in order to remove any traces of uncured resin. Alphaseal should be applied over conventional self-cured dental bases in case of direct pulp exposure.

INSTRUCTIONS FOR USE

1. Prepare the cavity in the conventional manner.
2. Rinse the cavity and dry.
3. Dispense Alphaseal onto a pad. After dispensing, slightly pull back the syringe plunger to reverse flow. **CAUTION:** shield the dispensed material from bright light before use.
4. Apply the material with a ballpoint applicator to provide a first coating of no more than 0.3-0.5mm thick. Remove any material inadvertently applied over enamel margins.
5. Cure for 20 seconds (30 seconds if the light is lower intensity), keeping the tip as close to the cured surface as possible.
6. Apply a second layer, if indicated. The second layer can be up to 1.5mm thick. Cure for 30 seconds.

STORAGE

Store at temperatures below 82°F/28°C. **DO NOT FREEZE.**

WARNING: Avoid prolonged contact of uncured Alphaseal with skin or gingivae.

Wash thoroughly if such contact occurs.

**Dent Zar, Inc.
19643 Trull Brook Dr.
Tarzana, CA 91356
U.S.A.**

Head Office

Toll Free: 800-444-1241, Outside U.S. 323-465-3621

Fax: 323-465-2994, E-mail: dentzar@yahoo.com., Web Site: www.dentzar.com.